



PRODUCT INFORMATION

Item Number	FN100
Product Colour	Blue
Size Range	S-XXL
Fingertip Thickness	0.15 mm
Palm Thickness	0.10 mm
Typical Weight	6.60 g
Shelf Life	5 Years

FEATURES

Chemical Resistant: Achieves type B chemical resistance according to EN ISO 374-1, for protection against at least 3 different chemicals

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria & fungi

Sensitivity: Latex Free; No latex proteins - eliminates protein sensitisation

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Powder Content: Powder Free, reducing the risk of dust contamination

Grip: Textured fingertips provide an improved grip and optimise surface performance in both wet and dry conditions

SECTORS

Automotive • Engineering and Manufacturing • Food Manufacturing and Preparation • Pharmaceutical and Laboratory



STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact



PPE CAT III



JKPT



VIRUS



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



PACKING

100 gloves per dispenser
10 dispensers per case
50 cases per pallet

OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5010699506955
Medium	5010699506962
Large	5010699506979
XL	5010699506986
XXL	5010699506993

Size	Outer Packaging Barcode
Small	05010699805218
Medium	05010699805225
Large	05010699805232
XL	05010699805249
XXL	05010699811950

Inner Packaging	
Inner dimensions (LxWxH)	240 x 130 x 70mm

Outer Packaging	
Outer Dimensions (LxWxH)	365 x 270 x 250mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555