

BIOPLUS

PRODUCT INFORMATION

Item Number	GL51
Product Colour	Cream
Size Range	S-XXL
Fingertip Thickness	0.20 mm
Palm Thickness	0.18 mm
Typical Weight	9.50 g
Shelf Life	5 Years

FEATURES

Fit: Ambidextrous - can be worn on either hand

Cuff: Beaded cuff to aid donning

Powder Content: Powder free, reducing the risk of dust contamination

Grip: Textured finish for improved grip

SECTORS

Medical, Vet, and Dental

STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425



VIRUS



STANDARDS



CAUTION

Contains latex which may cause allergic reactions, including anaphylactic response.
Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

BIOPLUS**PACKING**

100 gloves per bag
10 bags per case
32 cases per pallet

Size Outer Packaging Barcode

Small	5060040795538
--------------	----------------------

Medium	5060040795521
---------------	----------------------

Large	5060040795514
--------------	----------------------

XL	5060040795545
-----------	----------------------

XXL	5060040795552
------------	----------------------

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555