

PRODUCT INFORMATION

| | |
|---------------------|---------|
| Item Number | GL890 |
| Product Colour | Blue |
| Size Range | XS-XXL |
| Fingertip Thickness | 0.11 mm |
| Palm Thickness | 0.08 mm |
| Typical Weight | 5.00 g |
| Shelf Life | 5 Years |

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Sensitivity: Latex Free; No latex proteins - eliminates protein sensitisation

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Material Composition: Manufactured from nitrile butadiene polymer

Powder Content: Powder Free, reducing the risk of dust contamination

Grip: Textured finish to the fingers to enhance grip and handling

SECTORS

Automotive • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet, and Dental



STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact



MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



EN ISO 374-1:
2016 +A1:
2018 Type C



EN ISO 374-5:
2016



VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

PACKING

100 gloves per dispenser
10 dispensers per case
50 cases per pallet

INNER PACKAGING



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

OUTER PACKAGING



| Size | Inner Packaging Barcode |
|--------|-------------------------|
| XS | 5010699528322 |
| Small | 5010699528339 |
| Medium | 5010699528346 |
| Large | 5010699528353 |
| XL | 5010699528360 |
| XXL | 5010699536570 |

| Inner Packaging | |
|--------------------------|-------------------|
| Inner dimensions (LxWxH) | 235 x 117 x 78 mm |

| Size | Outer Packaging Barcode |
|--------|-------------------------|
| XS | 5010699839336 |
| Small | 5010699839343 |
| Medium | 5010699839350 |
| Large | 5010699839367 |
| XL | 5010699839374 |
| XXL | 5010699849809 |

| Outer Packaging | |
|--------------------------|-------------------|
| Outer Dimensions (LxWxH) | 418 x 258 x 251mm |

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555