

PRODUCT INFORMATION

Item Number	GL891
Product Colour	Blue
Size Range	S-XL
Fingertip Thickness	0.15 mm
Palm Thickness	0.10 mm
Typical Weight	7.80 g
Shelf Life	5 Years

FEATURES

Chemical Resistant: Achieves type B chemical resistance according to EN ISO 374-1, for protection against at least 3 different chemicals

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Sensitivity: Latex Free; No latex proteins - eliminates protein sensitisation

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Grip: Textured fingertips provide an improved grip and optimise surface performance in both wet and dry conditions

Cuff: Long cuff - 30cm glove for extra protection

Material Composition: Manufactured from nitrile butadiene polymer

Powder Content: Powder Free, reducing the risk of dust contamination

Suitability: Suitable for use with cytotoxic/chemotherapy drugs

SECTORS

Automotive • Engineering and Manufacturing • Food
Manufacturing and Preparation • Medical, Vet, and Dental



STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to ASTM D6978-05



MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



JKPT



VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

PACKING

100 gloves per dispenser
10 dispensers per case
45 cases per pallet

INNER PACKAGING



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5010699525901
Medium	5010699525918
Large	5010699525925
XL	5010699525932

Inner Packaging	
Inner dimensions (LxWxH)	270 x 135 x 80mm

Size	Outer Packaging Barcode
Small	5010699834980
Medium	5010699834997
Large	5010699835000
XL	5010699532017

Outer Packaging	
Outer Dimensions (LxWxH)	418 x 293 x 271mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555