

NITREX

PRODUCT INFORMATION

Item Number	GN16
Product Colour	Blue
Size Range	XS-XL
Fingertip Thickness	0.07 mm
Palm Thickness	0.05 mm
Typical Weight	3.20 g
Shelf Life	3 Years

FEATURES

Sensitivity: Accelerator free, Removing the risk of contact dermatitis caused by accelerator chemicals
Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Grip: Textured fingertips provide an improved grip and optimise surface performance in both wet and dry conditions

Powder Content: Powder Free, reducing the risk of dust contamination

Suitability: Suitable for use with cytotoxic/chemotherapy drugs

SECTORS

Automotive • Construction and Utility • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet, and Dental • Oil, Gas, and Renewables

STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified



STANDARDS

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact




MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



UK
CA



EN ISO 374-1:
2016 Type B



EN ISO 374-5:
2016



AQL
1.5

KPT

VIRUS



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

NITREX

PACKING

100 gloves per dispenser
10 dispensers per case
80 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



Size	Inner Packaging Barcode
XS	5060040797686
Small	5060040797693
Medium	5060040797709
Large	5060040797716
XL	5060040797723

Size	Outer Packaging Barcode
XS	5060040797631
Small	5060040797648
Medium	5060040797655
Large	5060040797662
XL	5060040797679

Inner Packaging	
Inner dimensions (LxWxH)	220 x 125 x 45mm

Outer Packaging	
Outer Dimensions (LxWxH)	260 x 245 x 235mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555