



PRODUCT INFORMATION

Item Number	GN31
Product Colour	Natural
Size Range	S-XL
Fingertip Thickness	0.11 mm
Palm Thickness	0.10 mm
Typical Weight	5.50 g
Shelf Life	5 Years

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN374-1, for chemical protection against at least 1 chemical

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Sensitivity: Chlorinated to reduce latex protein content

Material Composition: Flexible material for improved feel and fit reducing fatigue during extended wear
Manufactured from natural rubber latex

Powder Content: Powder free, reducing the risk of dust contamination

Grip: Textured Finish to the palm to improve grip and handling

SECTORS

Hair and Beauty • Janitorial • Medical, Vet, and Dental

STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)



STANDARDS

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4



MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



EN ISO 374-1:
2016 + A1:
2018 Type C



EN ISO 374-5:
2016



AQL
1.5



2777

VIRUS



CAUTION

Contains latex which may cause allergic reactions, including anaphylactic response.

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



PACKING

100 gloves per dispenser
10 dispensers per case
48 cases per pallet

OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951023548
Medium	5024951023555
Large	5024951023562
XL	5024951023579

Inner Packaging	
Inner dimensions (LxWxH)	225 x 125 x 75mm

Size	Outer Packaging Barcode
Small	5024951023593
Medium	5024951023609
Large	5024951023616
XL	5024951023623

Outer Packaging	
Outer Dimensions (LxWxH)	385 x 258 x 235mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555