

## NITREX

### PRODUCT INFORMATION

Item Number	GN34-MPL
Product Colour	White
Size Range	S-XL
Fingertip Thickness	0.21 mm
Palm Thickness	0.18 mm
Typical Weight	11.40 g
Shelf Life	5 Years

### FEATURES

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Powder Content:** Powder Free, reducing the risk of dust contamination

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Grip:** Textured finish for improved grip

### SECTORS

Pharmaceutical and Laboratory

### STANDARDS

AQL 1.5

Cleanroom Compatible - Class 100/ ISO 5

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested to European Standards EN420, EN ISO 374-1, EN 374-4 & EN ISO 374-5



VIRUS



### STANDARDS



### CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

### PACKING

100 gloves per bag  
10 bags per case  
24 cases per pallet

### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

## NITREX

Size	Inner Packaging Barcode
Small	5060040798911
Medium	5060040798928
Large	5060040798935
XL	5060040798942

Size	Outer Packaging Barcode
Small	5060040795613
Medium	5060040795620
Large	5060040795637
XL	5060040795644

Outer Packaging	
Outer Dimensions (LxWxH)	430 X 285 X 320 mm

**PLEASE NOTE:**

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555