

NITREX

PRODUCT INFORMATION

Item Number	GN51-MPL
Product Colour	White
Size Range	S-XL
Fingertip Thickness	0.18 mm
Palm Thickness	0.15 mm
Typical Weight	7.80 g
Shelf Life	5 Years

FEATURES

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Powder Content: Powder Free, reducing the risk of dust contamination

Grip: Textured finish for improved grip

SECTORS

Pharmaceutical and Laboratory

STANDARDS

AQL 1.5

Cleanroom Compatible - Class 100/ ISO 5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested to European Standards EN420, EN ISO 374-1, EN 374-4 & EN ISO 374-5



STANDARDS



VIRUS



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

100 gloves per bag
10 bags per case
36 cases per pallet

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

NITREX**DISPOSAL INSTRUCTIONS**

Dispose of the product responsibly in accordance with local regulations.

Size	Inner Packaging Barcode
Small	5060040796528
Medium	5060040796535
Large	5060040796542
XL	5060040796559

Size	Outer Packaging Barcode
Small	5060040796573
Medium	5060040796580
Large	5060040796597
XL	5060040796603

Outer Packaging	
Outer Dimensions (LxWxH)	315 X 305 X 330mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555