



#### PRODUCT INFORMATION

Item Number	GN830
Product Colour	Indigo
Size Range	S-XL
Fingertip Thickness	0.08 mm
Palm Thickness	0.06 mm
Typical Weight	5.00 g
Shelf Life	3 Years

#### FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Cuff:** Long cuff - 29cm glove for extra protection

**Sensitivity:** Accelerator free, Removing the risk of contact dermatitis caused by accelerator chemicals  
Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Material Composition:** Manufactured from nitrile butadiene polymer

**Powder Content:** Powder Free, reducing the risk of dust contamination

**Suitability:** Suitable for use with cytotoxic/chemotherapy drugs

**Grip:** Textured finish to the fingers to enhance grip and handling

#### SECTORS

Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet, and Dental



#### STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

**Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4**



MDR - Class I  
PPE CAT III



CLASS 1  
MEDICAL  
DEVICE



EN ISO 374-1:  
2016 Type C



EN ISO 374-5:  
2016



AQL  
1.5

2797

VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### INNER PACKAGING



Size	Inner Packaging Barcode
Small	5024951017073
Medium	5024951017080
Large	5024951017087
XL	5024951017103

Inner Packaging	
Inner dimensions (LxWxH)	260 x 126 x 70mm

#### PACKING

100 gloves per dispenser  
10 dispensers per case  
40 cases per pallet

#### OUTER PACKAGING



Size	Outer Packaging Barcode
Small	5024951017110
Medium	5024951017127
Large	5024951017134
XL	5024951017141

Outer Packaging	
Outer Dimensions (LxWxH)	365 x 260 x 270mm

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555