



#### PRODUCT INFORMATION

Item Number	GN83
Product Colour	Blue
Size Range	S-XL
Fingertip Thickness	0.12 mm
Palm Thickness	0.10 mm
Typical Weight	4.80 g
Shelf Life	3 Years

#### FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Material Composition:** Flexible material for improved feel and fit reducing fatigue during extended wear

**Grip:** Textured finish to the fingers to enhance grip and handling

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Material Composition:** Manufactured from nitrile butadiene polymer

**Powder Content:** Powder Free, reducing the risk of dust contamination

#### SECTORS

Automotive • Engineering and Manufacturing • Janitorial • Medical, Vet, and Dental



#### STANDARDS

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5



MDR - Class I  
PPE CAT III



CLASS 1  
MEDICAL  
DEVICE



EN ISO 374-1:  
2016+A1:  
2018 Type C



EN ISO 374-5:  
2016



AQL  
1.5



2777

VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

#### CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

#### PACKING

100 gloves per dispenser  
10 dispensers per case  
84 cases per pallet

#### INNER PACKAGING



#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951005025
Medium	5024951005032
Large	5024951005049
XL	5024951005056

Inner Packaging	
Inner dimensions (LxWxH)	230 x 120 x 63mm

Size	Outer Packaging Barcode
Small	5024951005063
Medium	5024951005070
Large	5024951005087
XL	5024951005094

Outer Packaging	
Outer Dimensions (LxWxH)	330 x 248 x 240mm

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555