



PRODUCT INFORMATION

Item Number	GN92
Product Colour	White
Size Range	XS-XL
Fingertip Thickness	0.09 mm
Palm Thickness	0.07 mm
Typical Weight	3.50 g
Shelf Life	5 Years

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Suitability: Suitable for use with cytotoxic/chemotherapy drugs

Material Composition: Flexible material for improved feel and fit reducing fatigue during extended wear
Manufactured from nitrile butadiene polymer

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Powder Content: Powder Free, reducing the risk of dust contamination

SECTORS

Automotive • Engineering and Manufacturing • Hair and Beauty • Janitorial • Medical, Vet, and Dental



STANDARDS

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5



MDR - Class I
PPE CAT III

2777



CLASS 1
MEDICAL
DEVICE



VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

PACKING

200 gloves per dispenser
10 dispensers per case
60 cases per pallet

INNER PACKAGING



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

OUTER PACKAGING



Size	Inner Packaging Barcode
XS	5024951012948
Small	5024951012955
Medium	5024951012962
Large	5024951012979
XL	5024951012986

Inner Packaging	
Inner dimensions (LxWxH)	220 x 120 x 82mm

Size	Outer Packaging Barcode
XS	5024951012993
Small	5024951013006
Medium	5024951013013
Large	5024951013020
XL	5024951013037

Outer Packaging	
Outer Dimensions (LxWxH)	425 x 250 x 230mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555