



## PRODUCT INFORMATION

Item Number	GN99
Product Colour	Blue
Size Range	S-XL
Fingertip Thickness	0.11 mm
Palm Thickness	0.08 mm
Typical Weight	5.40 g
Shelf Life	3 Years

## FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Strength:** Superior strength nitrile formulation provides good resistance to tears and punctures

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Material Composition:** Manufactured from nitrile butadiene polymer

**Powder Content:** Powder Free, reducing the risk of dust contamination

**Suitability:** Suitable for use with cytotoxic/chemotherapy drugs

## SECTORS

Automotive • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Medical, Vet, and Dental



## STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact



MDR - Class I  
PPE CAT III



CLASS 1  
MEDICAL  
DEVICE



EN ISO 374-1:  
2016 Type C



EN ISO 374-5:  
2016



AQL  
1.5

2777

VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

#### CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

#### PACKING

100 gloves per dispenser  
10 dispensers per case  
72 cases per pallet

#### INNER PACKAGING



Size	Inner Packaging Barcode
Small	5024951014690
Medium	5024951014706
Large	5024951014713
XL	5024951014720

#### Inner Packaging

Inner dimensions (LxWxH)	237 x 117 x 78mm
--------------------------	------------------

#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### OUTER PACKAGING



Size	Outer Packaging Barcode
Small	5024951014744
Medium	5024951014751
Large	5024951014768
XL	5024951014775

#### Outer Packaging

Outer Dimensions (LxWxH)	400 x 240 x 247mm
--------------------------	-------------------

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555