



PRODUCT INFORMATION

Item Number	MFNP100
Product Colour	Purple
Size Range	XS - XL
Fingertip Thickness	0.08 mm
Palm Thickness	0.06 mm
Typical Weight	4.20 g
Shelf Life	3 Years

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Fit: Unique shape delivers a glove with improved fit

Sensitivity: Accelerator free, Removing the risk of contact dermatitis caused by accelerator chemicals
Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Powder Content: Powder Free, reducing the risk of dust contamination

Suitability: Suitable for use with cytotoxic/chemotherapy drugs

SECTORS

Automotive • Construction and Utility • Engineering and Manufacturing • Food Manufacturing and Preparation • Janitorial • Oil, Gas, and Renewables



STANDARDS

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to European standards EN420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5



MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



EN ISO 374-1:
2016+A1:
2018 Type C



EN ISO 374-5:
2016



AQL
1.5

2777

VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



PACKING

100 gloves per dispenser
10 dispensers per case
60 cases per pallet

OUTER PACKAGING



Size	Inner Packaging Barcode
XS	5010699508973
Small	5010699508980
Medium	5010699508997
Large	5010699509000
XL	5010699509017

Size	Outer Packaging Barcode
XS	05010699816115
Small	05010699816122
Medium	05010699816139
Large	05010699816146
XL	05010699816153

Inner Packaging	
Inner dimensions (LxWxH)	230 x 120 x 68mm

Outer Packaging	
Outer Dimensions (LxWxH)	352 x 240 x 248mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555