



PRODUCT INFORMATION

Item Number	SD91E
Product Colour	Cool Blue
Size Range	S-XL
Palm Thickness	0.06 mm
Typical Weight	3.20 g
Shelf Life	5 Years

FEATURES

SafeDon: Packed in a Cuff 1st™ dispenser for reduced cross-contamination
Dispenses one glove at a time, eliminating glove wastage

Wall Mountable: Dispensers can be wall mounted with compatible glove dispenser holders

Grip: Textured finish to the fingertips, to enhance grip and tactile sensitivity

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Powder Content: Powder Free, reducing the risk of dust contamination

Chemical Resistant: Achieves type B chemical resistance according to EN ISO 374-1, for protection against at least 3 different chemicals

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

SECTORS

Automotive • Engineering and Manufacturing • Food Manufacturing and Preparation • Hair and Beauty • Janitorial • Medical, Vet, and Dental • Pharmaceutical and Laboratory



STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested to ASTM D6978-05

Tested to European Standard EN374-2

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact



MDR - Class I
PPE CAT III



CLASS 1
MEDICAL
DEVICE



KPT



VIRUS



2777



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

PACKING

200 gloves per dispenser
10 dispensers per case
72 cases per pallet

INNER PACKAGING



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951030065
Medium	5024951030072
Large	5024951030089
XL	5024951030096

Inner Packaging	
Inner dimensions (LxWxH)	128 x 120 x 110mm (S, M, L) / 128 x 120 x 130mm (XL)

Size	Outer Packaging Barcode
Small	5024951030027
Medium	5024951030034
Large	5024951030041
XL	5024951030058

Outer Packaging	
Outer Dimensions (LxWxH)	625 x 269 x 125mm (S, M, L) / 624 x 296 x 145mm (XL)

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555