



PRODUCT INFORMATION

Item Number	GN560
Product Colour	Pink
Size Range	S-L
Fingertip Thickness	0.09 mm
Palm Thickness	0.07 mm
Typical Weight	3.50 g
Shelf Life	5 Years

FEATURES

Chemical Resistant: Achieves type B chemical resistance according to EN ISO 374-1, for protection against at least 3 different chemicals

EN ISO 374-5 Protection: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Comfort: Nitrile technology provides an exceptionally soft and pliable material

Grip: Textured finish to the fingertips, to enhance grip and tactile sensitivity

Food Safe: Suitable for food contact; all food types. Please refer to our online Food Compatibility Chart for full details

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

Powder Content: Powder Free, reducing the risk of dust contamination

Silicone Free: Silicone residues from gloves can interfere with painting, coating, or adhesive bonding processes, potentially causing defects. These silicone-free gloves help minimise that risk, ensuring clean application surfaces.

SECTORS

Engineering and Manufacturing • Food Manufacturing and Preparation • Hair and Beauty • Janitorial • Medical, Vet, and Dental • Pharmaceutical and Laboratory • Retail and Household



STANDARDS

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

Tested in compliance with Commission Regulation (EU) No 10/2011 for Food contact

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4



MDR - Class I
PPE CAT III

2777



CLASS 1
MEDICAL
DEVICE



EN ISO 374-1:
2016+A1:
2018 Type B

KPT



EN ISO 374-5:
2016

VIRUS



UK
CA



AQL
1.5



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

100 gloves per dispenser
10 dispensers per case
90 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951016472
Medium	5024951016489
Large	5024951016496

Size	Outer Packaging Barcode
Small	5024951016502
Medium	5024951016519
Large	5024951016526

Inner Packaging	
Inner dimensions (LxWxH)	220 x 120 x 63 mm
Inner Packaging Weight	53g

Outer Packaging	
Outer Dimensions (LxWxH)	315 x 248 x 230 mm
Outer Packaging Weight	400g

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555