

PRODUCT INFORMATION

| | |
|---------------------|---------|
| Item Number | FT130 |
| Product Colour | Cream |
| Size Range | S-XL |
| Fingertip Thickness | 0.08 mm |
| Palm Thickness | 0.08 mm |
| Typical Weight | 7.10 g |
| Shelf Life | 5 Years |

FEATURES

Sensitivity: Latex free reducing the risk of hypersensitivity and allergic reactions

Powder Content: Powder free, reducing the risk of dust contamination

Comfort: A thin synthetic glove with an excellent fit that feels like a second skin
The soft vinyl formulation provides greater comfort and is more flexible and dextrous to wear

SECTORS

Automotive • Construction and Utility • Engineering and Manufacturing • Medical, Vet, and Dental

STANDARDS

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

100 gloves per dispenser
10 dispensers per case
40 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

INNER PACKAGING



| Size | Inner Packaging Barcode |
|--------|-------------------------|
| Small | 5010699507808 |
| Medium | 5010699507815 |
| Large | 5010699507822 |
| XL | 5010699507839 |

| Inner Packaging | |
|--------------------------|------------------|
| Inner dimensions (LxWxH) | 265 x 135 x 75mm |

OUTER PACKAGING



| Size | Outer Packaging Barcode |
|--------|-------------------------|
| Small | 05010699813947 |
| Medium | 05010699813954 |
| Large | 05010699813961 |
| XL | 05010699813978 |

| Outer Packaging | |
|--------------------------|-------------------|
| Outer Dimensions (LxWxH) | 350 x 278 x 283mm |

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555