



#### PRODUCT INFORMATION

Item Number	GD76
Product Colour	Clear
Size Range	M-L
Fingertip Thickness	0.05 mm
Palm Thickness	0.05 mm
Typical Weight	1.50 g
Shelf Life	5 Years

#### FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria & fungi

**Grip:** Embossed for improved grip

**Comfort:** Flexible material for improved feel and fit reducing fatigue during extended wear


**Material Composition:** Manufactured from Polyethylene

**Sensitivity:** Latex Free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Powder Content:** Powder free, reducing the risk of dust contamination

#### SECTORS

Hair and Beauty • Janitorial • Medical, Vet, and Dental



#### STANDARDS

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested to European Standards EN420, EN ISO 374-1, EN 374-4 & EN ISO 374-5

AQL 1.5



PPE CAT III



2777



#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



#### PACKING

250 gloves per dispenser  
10 dispensers per case  
96 cases per pallet

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### INNER PACKAGING

#### OUTER PACKAGING



Size	Inner Packaging Barcode
Medium	5024951013341
Large	5024951013358
XL	5024951013365

Size	Outer Packaging Barcode
Medium	5024951013372
Large	5024951013389
XL	5024951013396

Inner Packaging	
Inner dimensions (LxWxH)	244 x 120 x 55mm

Outer Packaging	
Outer Dimensions (LxWxH)	290 x 250 x 254mm

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555