

## PRODUCT INFORMATION

Item Number	GL621
Product Colour	Clear
Size Range	S-XL
Fingertip Thickness	0.06 mm
Palm Thickness	0.08 mm
Typical Weight	5.00 g
Shelf Life	5 Years

## FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN374-5, for protection against bacteria, fungi and virus

**Comfort:** A thin synthetic glove with an excellent fit that feels like a second skin

**Fit:** Ambidextrous. Can be worn on either hand

**Sensitivity:** Latex free reducing the risk of hypersensitivity and allergic reactions

**Powder Content:** Powder free, reducing the risk of dust contamination

## SECTORS

Janitorial • Medical, Vet, and Dental • Pharmaceutical and Laboratory

## STANDARDS

AQL 1.5

EU Certified

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

Tested to European Standards EN420 & EN ISO 374-1 & EN ISO 374-5



## STANDARDS

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425



MDR - Class I  
PPE CAT III



2777

VIRUS



## STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

## PACKING

100 gloves per dispenser  
10 dispensers per case  
84 cases per pallet

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### INNER PACKAGING



#### OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5010699528285
Medium	5010699528292
Large	5010699528308
XL	5010699528315

Inner Packaging	
Inner dimensions (LxWxH)	230 x 120 x 63mm

Size	Outer Packaging Barcode
Small	5010699839299
Medium	5010699839305
Large	5010699839312
XL	5010699839329

Outer Packaging	
Outer Dimensions (LxWxH)	333 x 248 x 240mm

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555