

## NITREX

### PRODUCT INFORMATION

Item Number	GN37
Product Colour	White
Size Range	S-XL
Fingertip Thickness	0.19 mm
Palm Thickness	0.17 mm
Typical Weight	11.40 g
Shelf Life	5 Years

### FEATURES

**Sterile:** Sterilised by irradiation method

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions

**Finish:** Soft smooth formulation for improved feel and fit reducing fatigue during extended wear

**Material Composition:** Manufactured from nitrile butadiene

**Powder Content:** Powder free, reducing the risk of dust contamination

**Grip:** Textured finish for improved grip

**Fit:** Hand specific to reduce hand fatigue

### SECTORS

Pharmaceutical and Laboratory

### STANDARDS

AQL 1.5

Cleanroom Compatible - Class 100/ ISO 5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested to European Standards EN420:2003+A1:2009, EN ISO 374-1:2016, EN ISO 374-4:2019 & EN ISO 374-5:2016



### STANDARDS



VIRUS



### CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

## NITREX

### PACKING

1 pair per pouch  
10 pouches per polybag  
20 polybags per case  
30 cases per pallet

### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

### INNER PACKAGING



### OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5060040792599
Medium	5060040792582
Large	5060040792575
XL	5060040792605

Size	Outer Packaging Barcode
Small	5060040799628
Medium	5060040799604
Large	5060040799611
XL	5060040799642

### Outer Packaging

**Outer Dimensions (LxWxH)**      **337 x 280 x 365mm**

### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555