

NITREX

PRODUCT INFORMATION

Item Number	GN38
Product Colour	White
Size Range	6-9
Fingertip Thickness	0.19 mm
Palm Thickness	0.16 mm
Typical Weight	16.50 g
Shelf Life	3 Years

FEATURES

Sterile: Sterilised by irradiation method

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions

Fit: Hand specific to reduce hand fatigue

Material Composition: Manufactured from nitrile butadiene

Grip: Textured finish for improved grip

SECTORS

Pharmaceutical and Laboratory

STANDARDS

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Tested to European Standards EN420 & EN ISO 374-1 & EN ISO 374-5

Cleanroom Compatible - Class 100/ ISO 5



VIRUS



STANDARDS



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

1 pair per pouch
25 Pairs per bag
4 Bags per case
35 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

NITREX

INNER PACKAGING



OUTER PACKAGING



Size	Inner Packaging Barcode
6	5060040792759
6.5	5060040792766
7	5060040792773
7.5	5060040792780
8	5060040792797
8.5	5060040792803
9	5060040792810

Size	Outer Packaging Barcode
6	5060040792827
6.5	5060040792834
7	5060040792841
7.5	5060040792858
8	5060040792865
8.5	5060040792872
9	5060040792889

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555