



## PRODUCT INFORMATION

Item Number	GN63
Product Colour	Cream
Size Range	S-XL
Fingertip Thickness	0.09 mm
Palm Thickness	0.07 mm
Typical Weight	4.50 g
Shelf Life	5 Years

## FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Comfort:** Soft smooth formulation for improved feel and fit reducing fatigue during extended wear

**Sensitivity:** Latex Free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Material Composition:** Manufactured from synthetic emulsion PVC polymer resin

**Powder Content:** Powder free, reducing the risk of dust contamination

## SECTORS

Hair and Beauty • Janitorial • Medical, Vet, and Dental

## STANDARDS

AQL 1.5

This product is classified as a Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified



## STANDARDS

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4


Tested to European standards EN420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5



MDR - Class I  
PPE CAT III



2777



VIRUS



## STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



#### PACKING

100 gloves per dispenser  
10 dispensers per case  
84 cases per pallet

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### INNER PACKAGING



#### OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951023838
Medium	5024951023845
Large	5024951023852
XL	5024951023869

Size	Outer Packaging Barcode
Small	5024951023876
Medium	5024951023883
Large	5024951023890
XL	5024951023906

Inner Packaging	
Inner dimensions (LxWxH)	230 x 120 x 63mm

Outer Packaging	
Outer Dimensions (LxWxH)	330 x 248 x 240mm

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555