



#### PRODUCT INFORMATION

Item Number	GN67
Product Colour	Green
Size Range	S - L
Fingertip Thickness	0.10 mm
Palm Thickness	0.08 mm
Typical Weight	6.10 g
Shelf Life	5 Years

#### FEATURES

**Coating:** Aloe Vera coating reacts to body heat, moisturising and hydrating skin to help keep hands soft and nourished

**Comfort:** Soft smooth formulation for improved feel and fit reducing fatigue during extended wear

**Sensitivity:** Latex Free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

**Material Composition:** Manufactured from synthetic emulsion PVC polymer resin

**Powder Content:** Powder free, reducing the risk of dust contamination

#### SECTORS

Hair and Beauty • Janitorial • Medical, Vet, and Dental

#### STANDARDS

AQL 1.5

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4



#### STANDARDS



#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

#### PACKING

100 gloves per dispenser  
10 dispensers per case  
84 cases per pallet

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

**INNER PACKAGING**



Size	Inner Packaging Barcode
Small	5024951006008
Medium	5024951006015
Large	5024951006022

Inner Packaging	
Inner dimensions (LxWxH)	230 x 120 x 63mm

**OUTER PACKAGING**



Size	Outer Packaging Barcode
Small	5024951006039
Medium	5024951006046
Large	5024951006053

Outer Packaging	
Outer Dimensions (LxWxH)	330 x 248 x 240mm

**PLEASE NOTE:**

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555