



## PRODUCT INFORMATION

Item Number	GN70
Product Colour	Blue
Size Range	S-XL
Fingertip Thickness	0.06 mm
Palm Thickness	0.05 mm
Typical Weight	5.40 g
Shelf Life	5 Years

## FEATURES

**Chemical Resistant:** Achieves type C chemical resistance according to EN374-1, for chemical protection against at least 1 chemical

**EN ISO 374-5 Protection:** Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

**Material Composition:** Hybrid - Latest technology combining strength of two materials

**Sensitivity:** Accelerator free, removing the risk of contact dermatitis caused by accelerator chemicals  
Latex free, reducing the risk of hypersensitivity and allergic reactions

**Material Composition:** Manufactured from flexible material for improved feel and fit reducing fatigue during extended wear  
Manufactured From a blend of Polyvinyl Chloride and Acrylonitrile butadiene

**Powder Content:** Powder Free - reduces risk of dust contamination

**Comfort:** The smooth finish offers good comfort and feel

## SECTORS

Automotive • Hair and Beauty • Medical, Vet, and Dental •  
Pharmaceutical and Laboratory



## STANDARDS

AQL 1.5

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with Medical Devices Regulation (EU) 2017/745 Class 1

Tested to European Standards EN420 & EN ISO 374-1 & EN ISO 374-5

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4



MDR - Class I  
PPE CAT III



2777

VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

#### INNER PACKAGING



#### PACKING

100 gloves per dispenser  
10 dispensers per case  
60 cases per pallet

#### OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951012122
Medium	5024951012146
Large	5024951012160
XL	5024951012184

Size	Outer Packaging Barcode
Small	5024951012139
Medium	5024951012153
Large	5024951012177
XL	5024951012191

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555