



#### PRODUCT INFORMATION

Item Number	GS641
Product Colour	Clear
Size Range	S-L
Fingertip Thickness	0.03 mm
Palm Thickness	0.03 mm
Typical Weight	2.30 g
Shelf Life	5 Years

#### FEATURES

**Sterile:** Sterilised by irradiation method

**Sensitivity:** Latex free, reducing the risk of hypersensitivity and allergic reactions

**Comfort:** Flexible Material for improved feel and fit reducing fatigue during extended wear

**Material Composition:** Manufactured from evathene

**Cuff:** Non Beaded to allow for easy over glove donning and removal

**Powder Content:** Powder free, reducing the risk of dust contamination

**Finish:** Smooth finish for perfect outer glove donning and sensitivity

#### SECTORS

Medical, Vet, and Dental • Pharmaceutical and Laboratory

#### STANDARDS

AQL 1.5

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

This product is classified as Class Is Medical Device in compliance with the requirements of Medical Directive 93/42/EEC on Medical Devices (MDD), Annex V



#### STANDARDS



STERILE R



2460

**UK  
CA**

#### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

#### PACKING

1 pair per pouch  
50 pairs per dispenser  
24 dispensers per case  
25 cases per pallet

#### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

## INNER PACKAGING



## OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951024446
Medium	5024951024453
Large	5024951024460

Size	Outer Packaging Barcode
Small	5024951024477
Medium	5024951024484
Large	5024951024491

Inner Packaging	
Inner dimensions (LxWxH)	115 x 100 x 175mm

Outer Packaging	
Outer Dimensions (LxWxH)	620 x 240 x 370mm

### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555