



PRODUCT INFORMATION

Item Number	GS690
Product Colour	Blue
Size Range	S-L
Fingertip Thickness	0.10 mm
Palm Thickness	0.08 mm
Typical Weight	4.20 g
Shelf Life	3 Years

FEATURES

Sterile: Sterilised by irradiation method

Strength: Offers better puncture tear and abrasion resistance than natural rubber or vinyl

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions

Comfort: Flexible Material for improved feel and fit reducing fatigue during extended wear

Material Composition: Manufactured from nitrile butadiene

Powder Content: Powder free, reducing the risk of dust contamination

Grip: Textured finish to the fingertips, to enhance grip and tactile sensitivity

SECTORS

Medical, Vet, and Dental • Pharmaceutical and Laboratory

STANDARDS

AQL 1.5

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

This product is classified as Class Is Medical Device in compliance with the requirements of Medical Directive 93/42/EEC on Medical Devices (MDD), Annex V



STANDARDS



STERILE R



**UK
CA**

CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

1 pair per pouch
50 pairs per dispenser
4 dispensers per case
56 cases per pallet

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



OUTER PACKAGING



Size	Inner Packaging Barcode
Small	5024951024248
Medium	5024951024255
Large	5024951024262

Inner Packaging	
Inner dimensions (LxWxH)	140 x 130 x 285mm

Size	Outer Packaging Barcode
Small	5024951024279
Medium	5024951024286
Large	5024951024293

Outer Packaging	
Outer Dimensions (LxWxH)	531 x 147 x 292mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555