

PRODUCT INFORMATION

Item Number	GSSYN
Product Colour	Natural
Size Range	5.5-9.0
Fingertip Thickness	0.23 mm
Palm Thickness	0.20 mm
Typical Weight	11.00 g
Shelf Life	3 Years

FEATURES

Sterile: Sterilised using electron beam > 25kGy

Fit: Inner Surface treated for ease of donning

Grip: Micro Roughened to provide excellent feel and instrument grip in both wet and dry conditions

Comfort: Manufactured on ergonomically designed formers from an advanced synthetic formulation, offers ultimate comfort and tactile sensitivity
Intrinsically soft with good elasticity.

Ergonomically designed with curved fingers for comfort and minimised hand fatigue
Flexible Material for improved feel and fit reducing fatigue during extended wear

Fit: Hand specific to reduce hand fatigue

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions

Powder Content: Powder free, reducing the risk of dust contamination

SECTORS

Medical, Vet, and Dental



STANDARDS

AQL 0.65

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

This product is classified as Class IIa Medical Device in compliance with the requirements of Medical Directive 93/42/EEC on Medical Devices (MDD)



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

PACKING

1 pair per pouch
50 pairs per dispenser
4 dispensers per case
36 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



Size	Inner Packaging Barcode
5.5	5024951024781
6	5024951024798
6.5	5024951024804
7	5024951024811
7.5	5024951024828
8	5024951024835
8.5	5024951024842
9	5024951024859

Size	Outer Packaging Barcode
5.5	5024951024866
6	5024951024873
6.5	5024951024880
7	5024951024897
7.5	5024951024903
8	5024951024910
8.5	5024951024927
9	5024951024934

Inner Packaging	
Inner dimensions (LxWxH)	270 x 140 x 250mm

Outer Packaging	
Outer Dimensions (LxWxH)	458 x 147 x 508mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555