



## PRODUCT INFORMATION

Item Number	INTW10
Product Colour	Natural
Size Range	5.5-9
Fingertip Thickness	0.23 mm
Palm Thickness	0.20 mm
Typical Weight	12.70 g
Shelf Life	5 Years

## FEATURES

**Sterile:** Sterilised by irradiation method

**Fit:** Hand specific to reduce hand fatigue

**Grip:** Micro Roughened to provide excellent feel and instrument grip in both wet and dry conditions

**Comfort:** Manufactured on ergonomically designed formers from an advanced synthetic formulation, offers ultimate comfort and tactile sensitivity

**Cuff:** Long beaded cuff

**Material Composition:** Manufactured from natural rubber latex, which can cause allergic reactions, including anaphylactic responses.

## SECTORS

Medical, Vet, and Dental

## STANDARDS

AQL 0.65

This product is classified as Class IIa Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No. 618 as amended) (UKCA)

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425



## STANDARDS

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5



PPE CAT III

STERILE R



KPT



VIRUS



## CAUTION

**Contains latex which may cause allergic reactions, including anaphylactic response.**

**Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.**

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

### INNER PACKAGING



Size	Inner Packaging Barcode
5.5	9555256608049
6	9555256608056
6.5	9555256608063
7	9555256608070
7.5	9555256608087
8	9555256608094
8.5	9555256608100
9	9555256608117

### PACKING

50 gloves per dispenser  
4 dispensers per case  
48 cases per pallet

Size	Outer Packaging Barcode
5.5	9555256607967
6	9555256607974
6.5	9555256607981
7	9555256607998
7.5	9555256608001
8	09555256608018
8.5	9555256608025
9	9555256608032

### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555