



PRODUCT INFORMATION

Item Number	INTW82E
Product Colour	Blue
Size Range	S-L
Fingertip Thickness	0.12 mm
Palm Thickness	0.10 mm
Typical Weight	6.16 g
Shelf Life	3 Years

FEATURES

Sterile: Sterilised by irradiation method

Suitability: Suitable for use with cytotoxic/chemotherapy drugs

Grip: Textured finish for improved grip

Cuff: Beaded extended cuff for extra protection and to aid donning

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions to natural rubber proteins

SECTORS

Medical, Vet, and Dental • Pharmaceutical and Laboratory

STANDARDS

AQL 1.5

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

This product is classified as Class Is Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No. 618 as amended) (UKCA)



STERILE R



2797



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

50 gloves per dispenser
6 dispensers per case
32 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



INNER PACKAGING



Size	Inner Packaging Barcode
Small	9555256603761
Medium	9555256603778
Large	955 5256603785

Inner Packaging	
Inner dimensions (LxWxH)	250x 150 x 135mm

Size	Outer Packaging Barcode
Small	09555256603815
Medium	09555256603822
Large	09555256603839

Outer Packaging	
Outer Dimensions (LxWxH)	460 x 280 x 260mm

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555