

## PRODUCT INFORMATION

Item Number	SDGKF
Product Colour	Green
Size Range	5.5-9.0
Fingertip Thickness	0.18 mm
Palm Thickness	0.17 mm
Typical Weight	10.30 g
Shelf Life	3 Years

## FEATURES

**Colour:** Green colour for clear detection of perforations when double gloving

**Grip:** Enhanced tactile grip provides excellent feel and instrument grip in both wet and dry conditions

**Fit:** Hand specific to reduce hand fatigue

**Material Composition:** Manufactured from natural rubber latex, which can cause allergic reactions, including anaphylactic responses.

**Powder Content:** Powder free, reducing the risk of dust contamination

**Sterile:** Sterilised using electron beam > 25kGy

**Finish:** Polymer Coated

## SECTORS

Medical, Vet, and Dental

## STANDARDS

AQL 0.65

This product is classified as Class IIa Medical Device in compliance with the requirements of Medical Directive 93/42/EEC on Medical Devices (MDD)

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425



## STANDARDS

**Compliant with European Standards EN455:** Medical Gloves for Single Use parts 1,2,3&4

Tested to European Standards EN420 & EN ISO 374-1 & EN ISO 374-5



KPT

VIRUS



## CAUTION

Contains latex which may cause allergic reactions, including anaphylactic response.

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

**STORAGE INSTRUCTIONS**

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

**PACKING**

1 Pair Per Pouch  
50 Pouches Per Dispenser  
4 Dispensers per Case

**DISPOSAL INSTRUCTIONS**

Dispose of the product responsibly in accordance with local regulations.

Size	Outer Packaging Barcode
5.5	19555076112716
6	19555076112723
6.5	19555076112730
7	19555076112747
7.5	19555076112754
8	19555076112761
8.5	19555076112778
9	19555076112785

**PLEASE NOTE:**

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555