

PRODUCT INFORMATION

| | |
|---------------------|---------|
| Item Number | SNA |
| Product Colour | Cream |
| Size Range | 5.5-9 |
| Fingertip Thickness | 0.20 mm |
| Palm Thickness | 0.18 mm |
| Typical Weight | 14.00 g |
| Shelf Life | 3 Years |

FEATURES

Sterile: Sterilised by irradiation method

Sensitivity: Accelerator free synthetic gloves removing the risk of contact dermatitis caused by these chemicals

Grip: Micro Roughened to provide excellent feel and instrument grip in both wet and dry conditions

Comfort: Manufactured on ergonomically designed formers from an advanced synthetic formulation, offers ultimate comfort and tactile sensitivity

Fit: Hand specific to reduce hand fatigue

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions

SECTORS

Medical, Vet, and Dental

STANDARDS

AQL 0.65

This product is classified as Class IIa Medical Device in compliance with the requirements of Medical Directive 93/42/EEC on Medical Devices (MDD)

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425



STANDARDS

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to European Standards EN420 & EN ISO 374-1 & EN ISO 374-5



PPE CAT III

STERILE R



CLASS IIa
MEDICAL
DEVICE



KMNPST



VIRUS



**UK
CA**

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

PACKING

1 pair per pouch
50 gloves per dispenser
4 dispensers per case
30 cases per pallet

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

INNER PACKAGING



| Size | Inner Packaging Barcode |
|------|-------------------------|
| 5.5 | 5060040798621 |
| 6 | 5060040798638 |
| 6.5 | 5060040798645 |
| 7 | 5060040798652 |
| 7.5 | 5060040798669 |
| 8 | 5060040798676 |
| 8.5 | 5060040798683 |
| 9 | 5060040798690 |

| Inner Packaging | |
|--------------------------|-------------------|
| Inner dimensions (LxWxH) | 290 x 153 x 225mm |

| Size | Outer Packaging Barcode |
|------|-------------------------|
| 5.5 | 5060040798539 |
| 6 | 5060040798546 |
| 6.5 | 5060040798553 |
| 7 | 5060040798560 |
| 7.5 | 5060040798577 |
| 8 | 5060040798584 |
| 8.5 | 5060040798591 |
| 9 | 5060040798614 |

| Outer Packaging | |
|--------------------------|-------------------|
| Outer Dimensions (LxWxH) | 595 x 315 x 235mm |

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555