



### PRODUCT INFORMATION

Item Number	GL2950
Product Colour	Blue
Size Range	One Size
Fingertip Thickness	0.02 mm
Palm Thickness	0.02 mm
Palm Width	180mm
Typical Weight	7.20 g
Shelf Life	5 Years

### FEATURES

**Length:** 900mm long

**Thickness:** 25µ thick

**Material Composition:** Flexible material for improved feel and fit reducing fatigue during extended wear

**Sensitivity:** Latex free reducing the risk of hypersensitivity and allergic reactions

**Material Composition:** Manufactured from Polyethylene

**Powder Content:** Powder free, reducing the risk of dust contamination

**Finish:** The smooth finish offers good comfort and feel

### SECTORS

Medical, Vet, and Dental

### STANDARDS

**CE Category:** Cat I Minimal Risk



### STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

### PACKING

100 gloves per dispenser  
20 dispensers per case  
50 cases per pallet

### DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



### OUTER PACKAGING



Size	Outer Packaging Barcode
Large	5055066587434

Inner Packaging	
Inner dimensions (LxWxH)	340 x 250 x 20mm

Outer Packaging	
Outer Dimensions (LxWxH)	470 x 260 x 215mm

#### PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555